

## Simplified Blood Lactate Test Meter

# Lactate Pro™ LT-1710

A blood lactate test meter in the size of a credit card, the "Lactate Pro" speedily measures lactate with only a small sample of blood. Realizing the C. V. 3%, "Lactate Pro" is proud of its high precision and performance, equivalent of large-scale testing instruments.



**Only a 5µL blood sample is required.**

**Speedy measurements in 60 sec.**

**High performance in a small size.**



- Required blood sample volume is 5µL.
- Measurements completed in 60 seconds.
- High data precision of C.V. 3%.
- The size of a credit card, the smallest in the world.
- Stores up to 20 measurement results and can average past results.

# Lactate Pro™ LT-1710

Simplified Blood Lactate Test Meter

**With only 5µL of blood and simple operation, "Lactate Pro" provides accurate blood lactate level in 60 seconds.**

"Lactate Pro" adopts the electrochemical measurement method utilizing enzymic reaction. The operation has been further simplified. As little as 5µL of blood is enough for a measurement. Furthermore, the wiping of excess blood and other troublesome

preparation are not required. Measurements are completed in a mere 60 seconds. The automatic aspiration of blood greatly reduces the risk of human error, and thus provides accurate blood lactate level without involving any manual operation.

■ Speedy measurements and simple operation



1 Insertion of the test strip.



2 Collection of blood.



3 Start of measurement.



4 Display of blood lactate level.

- 1 Insert the test strip into the strip inlet of the meter.
- 2 Before collecting 5µL of blood from your finger, wash and dry your hands.
- 3 Blood is automatically aspirated and measurement begins.
- 4 The measurement result is displayed in 60 seconds.


**High data precision of C.V. 3% satisfies the needs of clinical testing.**

With the high repeatability of C.V. 3%, "Lactate Pro" is able to be used not only in physical training but also useful for clinical testing. Thus, "Lactate Pro" provides reliably precise data equivalent to that of large-scale testing equipment and also maintains good correlation with conventional lactate test meters.

■ Repeatability test

n	Sample1	Sample2
1	2.1	11.2
2	2.3	11.4
3	2.3	11.4
4	2.1	11.4
5	2.2	11.3
6	2.2	11.7
7	2.3	11.4
8	2.2	11.4
9	2.2	11.9
10	2.2	12.0
11	2.1	11.4
12	2.2	12.0
13	2.1	12.1
14	2.2	11.8
15	2.2	11.8
mean	2.2	11.6
S.D.	0.07	0.39
C.V.	3.2	3.3

■ Correlation test



**A credit card size meter is very convenient for mobile use.**

Measuring 83.8 x 55.0 x 14.5 mm, the smallest of its kind in the world, "Lactate Pro" requires only a limited space. Further, weighing only 50g, it is easy to carry about. Whenever and wherever necessary, accurate blood lactate levels can be obtained.

**Memory function for the benefit of data management. Up to 20 measurement results can be stored.**

From the latest test result, up to 20 previous measurements can be recorded. Furthermore, the average of the past 20 measurement results can be displayed for the benefit of better data management and analysis.

## Specifications

### ● Meter

Test	Lactate in whole blood
Sample size	Approximately 5µL
Test strip	Exclusive Lactate Pro™ Test Strip
Measuring range	0.8~23.3mmol/L
Measuring time	60 seconds
Temperature compensation	Automatic compensation using a built-in sensor
Calibration	Automatically selects the appropriate calibration curve by using a calibration strip
Power source	Two 3-volt lithium batteries (DL or CR2032)

Battery life	Approximately 1000 tests
Operating environment	Temperature range 10°-40°C (50°-104°F), Humidity 20-80%RH
Dimensions	83.8(L) × 55.0(W) × 14.5(H) mm
Weight	50g

### ● Test Strip

Storage	Store at room temperature lower than 30°C
Expiry	1 year after production
Package unit	25 strips/box

Design and specifications may be subject to change without notice.



57 NISHI AKETA-CHO, HIGASHI-KUJO, MINAMI-KU,  
KYOTO 601-8045, JAPAN  
TEL +81-75-662-8967 FAX +81-75-662-8973

Distributed By:  
Vetlab Supply  
Palmetto Bay, FL 33157  
(800)330-1522 Fax (305)232-8421  
www.vetlab.com