


IRMA TRUPOINT®

Blood Analysis System



*Accurate, reliable results
where and when you need them.*

When managing the most severe and life threatening emergencies, speed, reliability and accuracy is everything. IRMA TRUpoint sets the new standard for rapid testing of blood gas, electrolytes and other time critical parameters in the veterinary clinic or in the field. This powerful, self-contained system can go wherever you need it for immediate results and more timely medical intervention.



ACCURATE

Obtain immediate results for improved efficiencies, diagnostic capabilities and procedural outcomes.

- Precise results within two minutes
- Lab-quality accuracy
- Comprehensive menu
- Custom correlation for site-specific analyte reference ranges
- Ability to establish species specific references ranges

RELIABLE

Trust, integrity and performance with simple, intuitive operation and advanced quality features.

- Single-use cartridge for maintenance-free convenience and reproducibility
- Self-contained calibrant in each cartridge
- Integrated, automatic electronic quality control for enhanced safety and improved compliance
- Exceptional service and support

EASY TO USE

Designed for ease of use, the IRMA TRUpoint can be used by anyone, anywhere.


- Complete portability with lightweight, easy grip, battery operation
- Intuitive, touch screen software guides operator and is easy to learn
- Room temperature cartridges are ready to test and go at anytime
- Cartridge based, closed system simplifies sample introduction





SINGLE-USE CARTRIDGE: ANOTHER TRUPOINT ADVANTAGE

An important feature of the IRMA TRUpoint system, is the single-use cartridges which are available in a variety of analyte configurations. Each cartridge contains miniaturized electrochemical measurement sensors and an onboard calibrant. With a wide operating test range, IRMA TRUpoint can be used in colder environments. Most cartridges do not require refrigerated storage and are ready to use anywhere, anytime.



ITC

For more than 36 years, ITC has been at the Point of Care. By focusing exclusively on point of care testing longer than anyone else, we've always made a difference wherever you need it. And with service and support that's second to none, we're by your side, so you can stay by theirs.

To learn more about IRMA TRUpoint and how ITC can help you improve patient care, please contact us at 1.800.631.5945 or visit www.itcmed.com/irma.shtml.

©2006 ITC. All rights reserved. IRMA TRUpoint is a trademark and IRMA is a registered trademark of ITC.

CARTRIDGE TEST MENU

	CC	BG	H3	H4	GL	CR	LA
Measured Values	pH	*	*				
	pCO ₂	*	*				
	pO ₂	*	*				
	Na ⁺	*		*	*	*	
	K ⁺	*		*	*	*	
	Cl ⁻				*	*	
	iCa	*		*			
	Glu					*	
	BUN				*		
	Hct	*		*	*		
Calculated Values	Creatinine					*	
	Lactate						*
	HCO ₃ ⁻	*	*				
	TCO ₂	*	*				
	BE _{ecf}	*	*				
	BE _b	*	*				
	O ₂ SAT	*	*				
	iCa(N)	*					
	tHb	*		*	*		



Distributed By:
Vetlab Supply
(800)330-1522
(305)232-8421 Fax
Visit us online at www.vetlab.com