

# MobilPro™

## Cordless Rechargeable Microscope

The MobilPro is a cordless microscope that goes wherever you need it! Powered by a rechargeable NiMH battery that delivers 50-70 hours of continuous use between charges, the MobilPro uses cool LED illumination instead of a traditional bulb to provide you with clear, sharp images. The unit can be fully recharged in just 8 hours and batteries can be cycled up to 500 times before needing replacement. The MobilPro comes complete with dust cover, battery, and recharging station.



- 4 Position Reverse Nosepiece with 4X,10X, 40X, and 100X Oil Brightfield Objectives.
- Cool LED Illumination. *(No bulbs to change!)*
- Long Life NiMH Battery Recharges in just 8 hours.
- 50-70 Hours of Use from a Single Charge.

