



# MobilPro™

## Cordless Rechargeable Microscope Goes Anywhere, Anytime!



- **4 Position Reverse Nosepiece with 4X, 10X, 40X, and 100X Oil Brightfield Objectives**
- **Cool LED Illumination (No bulbs to change!)**
- **NiMH Battery Recharges in just 8 hours**
- **50-70 Hours of use from a Single Charge**

*The MobilPro is a cordless binocular microscope that goes wherever you need it. Powered by a rechargeable NiMH battery that delivers 50-70 hours of continuous use between charges, the MobilPro uses cool LED illumination instead of a traditional bulb to provide you with clear, sharp images. The unit can be fully recharged in just 8 hours and the batteries can be cycled up to 500 times before needing replacement. The MobilPro comes complete with dust cover, battery, and recharging station.*

**Introductory Special! Just 669.00**



*Distributed By:*

**Vetlab® Supply, Inc.**

*"The In-house Laboratory Specialist"*

18131 SW 98th Ct.

Miami, FL 33157

305-253-1848 · 800-330-1522

Fax 305-232 8421

Visit us online at [www.vetlab.com](http://www.vetlab.com)