

#### www.unicosci.com

# PowerSpin™ DX Centrifuge

## **Product Images**







**Short Description** 

POWERSPIN™ DX CENTRIFUGE VARIABLE SPEED 1000-3400 RPM VARIABLE SPEED 500-3500 RPM PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL 6 PLACE HORIZONTAL ROTOR 60 MIN DIGITAL TIMER

## **Description**

POWER INPUT 110V

CAPACITY 6X15 ML TUBES

UNICO is proud to present the PowerSpin™ DX Centrifuge

With 3 rotor options this centrifuge is a great answer for every day laboratory tasks. Choose from 24 place fixed rotor for 2-15ml tubes (#C8724), 6 place horizontal rotor for 2-15ml tubes (#C8706) or 4 place horizontal rotor for up to 50ml tubes (#C8704)

The PowerSpin DX series features a large, spacious stainless steel lined interior, all digital push-button control with 3 preset operational modes for Blood, Urine and Fecal, variable speed 500-4000 RPM and zero RPM lid latching system. DX's are available in three main versions, fixed angle and horizontal spin rotors. All DX's feature maintenance-free brushless motors and suction cup rubber feet.

### Additional Information

	Header	Acces	ssories
Accessories	C8600-01	6 Plac	ce Horizontal Rotor, 2-15ml for DX Centrifuge
	C8600-02	4 Plac	ce Horizontal rotor, 50ml DX Centrifuge
	C8600-03M	24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)	
	C800-02	Regul	ar tube shield, pack of 2
	C800-03	Adap	ter for short tubes, pack of 10
	C800-04	Pedia	tric or microtube adapter, pack of 2
	C800-05	Tube	Lift Adapters, Pk/12
Specification	Speed (RPM	1)	500~4000 rpm
	Maximum RCF (g) 2508		
	Rotor		6 place horizontal
	Tube Size		2-15ml
	Motor		DC
	Power Inpu	t	110V
	Size		13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width
Size and Weight	Shipping Siz	ze '	19"X18"X15"

Weight

Shipping Weight 36 LBS

**31 LBS**