



# PowerSpin™ MX Horizontal Spin Centrifuge

UNICO® is proud to introduce our powerful 6-place Horizontal Spin centrifuge, the new PowerSpin™ MX! The PowerSpin™ MX spins 1000 – 3400 RPM, to insure a clean and complete separation of cellular components yielding superior serum or plasma specimen. Horizontal swing-out rotor provides high G-force and flat separation meniscus line.

The all digital push button front panel makes it easy to select or modify speed and run time even while wearing gloves. The PowerSpin™ MX also features a two line digital display that shows status, time, speed (RPM), and RCF (g).

Each PowerSpin™ Series centrifuge runs extremely quiet. The rubber shock absorbing suction cup feet hold the unit in place resulting in a quieter run with less vibration

Integrated True Zero RPM Safety System prevents lid from being open while the rotor is in motion.

Call us today for more information.



### Features and Specifications:

- 6 place swing-out horizontal spin rotor for 2 to 15 ml tubes
- Push button panel for convenience when wearing gloves
- Quiet running , only 50 decibels at 3400 RPM
- Stainless steel rotor chamber is accessible for easy cleaning
- Easy to read digital display shows real time RPM and RCF readings
- Remembers last setting for quick and easy repeat testing
- Variable Speed, 1000-3400 RPM
- Suction cup rubber feet eliminate vibration
- 60 minute digital timer with tone signal
- On/off switch located on rear panel of unit
- Detachable power cord
- Lid window for speed verification with strobe tachometer
- Zero RPM Safety Latching System to prevent lid from being opened while rotor is spinning
- Manual lock override in case of power outage
- Indicator light signals when the power is on
- Standard 115V/60Hz AC (220V/50Hz available on request)

### Ordering Information:

Model	Catalog #	Speed Range	Speed Monitor	Max. RCF	Capacity	Swing-out Rotor	Motor
<b>PowerSpin™ MX</b>	<b>8606</b>	1000 - 3400 RPM	Digital display	1825 G	6 x 2-15 ml	6 place	AC 115V

### Accessories and parts:

C8606-01	6 place swing out rotor	Centrifuge Size:	13 3/4" x 15" x 11"
C800-02	Regular tube shield, pack of 2	Shipping Size:	18 in x 19 in x 15.5"
C800-03	Adapters for 5 ml tubes, pack of 10	Net Weight:	31 lbs
C800-04	Pediatric or 2 ml tube adapters, pack of 2	Shipping Weight:	36 lbs
C800-10	Fuse 250A 4A, each	Warranty:	One year, parts and labor

**United Products & Instruments, Inc.**  
182 Ridge Road, Suite E, Dayton, New Jersey 08810  
© 2005 UNICO®

Distributed by:  
**Vetlab Supply**  
Palmetto Bay, FL 33157  
(800)330-1522 • (305)232-8421  
www.vetlab.com