

In-Office Screening for Urinary Tract Infection (UTI) is as Easy as 1-2-3!



Uricult® dip paddles are an easy and cost effective way to detect common uropathogens in only 18-24 hours. Just apply sample to the paddle, incubate, and compare results to the color chart provided.

It's just that simple!

Test now, not later.

Infection of the urinary tract (UTI) is one of the most common disorders of dogs and cats. While many infections respond to a single course of broad spectrum antibiotics, obtaining a urine culture before treatment is initiated is just better medicine. The Uricult® dip paddle provides veterinarians with a convenient and cost effective tool that can be utilized at the initial visit.

Why Uricult®?

Simplicity. The Uricult® paddle is formulated with two specialized agars, CLED and EMB. During incubation, the CLED agar undergoes a color change that's dependent upon the organism present. EMB is a Gram negative selective agar that is especially helpful for identifying E. coli. Presumptive identification of many common uropathogens can be made by comparing the growth present on the Uricult® paddle to the color interpretation chart provided. A colony count may also be performed to assess the severity of infection. Cultures that are positive can be sent to the reference laboratory for confirmation and antibiotic susceptibility testing.



**Uricult C.L.E.D./EMB
Product No.ODU1000**

Distributed by:

Vetlab Supply • Palmetto Bay, FL • 800.330.1522 • Fax 305.232.8421 • www.vetlab.com