

# Video and Digital Microscopy

**Vetlab Supply offers a variety of full-featured video cameras for many applications and to fit any budget. All units can be eyepiece mounted or affixed to a trinocular head for simultaneous ocular and video viewing.**



**Video Eyepiece Camera.** Economical basic camera that will project images on any video monitor or television equipped with standard RCA inputs. 640 X 480 resolution will produce very good images on smaller (<20") TV monitors. Images can also be output to VCR and recordable DVD equipment. Comes with built in 10X eyepiece and features automatic fine focus. Comes complete with power pack and 3 eyepiece adapters to fit most standard compound microscopes. **New lower price! 199.00**

**USB Eyepiece Camera Model SC3000** Budget-priced USB camera for capturing digital images to any Windows based computer. At 30 frames per second, USB 2.0 systems can capture near real time images at 3 megapixels. The SC3000 slides into most standard eyepiece tubes. Optional adapters are also available. Complete with Scopephoto Windows® capture software and USB cable. **435.00 (5.0 megapixel mod. SC5000 available 635.00)** \*Case design subject to change.



**Capture High Resolution Images with the 318CU Digital Camera.** The 3 mega-pixel 318CU is an easy to use quality camera for capturing high-resolution digital images and video from any microscope. Images are uploaded to computer via a high speed plug and play USB 2.0 connection. Included with the 318CU is Micrometrics' robust imaging software that lets you view real time full screen video while simultaneously capturing still images. Captured Images can be saved in any of 10 file formats including JPEG, TIFF, and BMP. Video is saved as an AVI file. The camera accepts most standard C-mount adapters. **899.00**



**Dual Output Moticam 480 Offers a Choice of Viewing Modes.** The Moticam 480 lets you view real time images on your TV monitor while simultaneously capturing digital images on your computer with just one click of the mouse. The Moticam's 1/3" CCD chip provides excellent quality 480 line video images that can be output in either standard or S-Video mode for even finer detail. Hi speed USB 2.0 delivers 640 X 480 pixel images. The Moticam 480 is priced at **879.00** and comes complete with adapters to fit almost any microscope, S-Video, USB, and RCA, cables, power supply, Windows® and Mac based



imaging software, and instructions. **\*The Moticam 480 may also be used alone for capturing and viewing macro images with an optional macro adapter and gooseneck stand.**

Need more information? Contact us at:

**Vetlab® Supply**

Palmetto Bay, FL 33157

(800) 330-1522 · (305) 253-1848 · (305) 232-8421

or visit us online at [www.vetlab.com](http://www.vetlab.com)

