


Bullseye Culture Systems

Features

- *Four large numbered identification sectors*
- *Full color wall charts for easy organism identification*
- *Complete and easy to follow instructions*
- *Large central area for antibiotic sensitivity test*
- *Sensitivity zone interpreter*
- *Designed for no open flames – disposable calibrated loops*
- *CO₂ atmosphere– easily achieved when necessary*

Benefits

- *Accurate and reliable results within 18-24 hours*
- *Low cost – Saves time and money*
- *No special training required – system can be utilized by current office personnel*
- *Choice of effective antibiotic made easy*
- *Completely disposable*
- *Toll free technical service available (MT, ASCP)*



Bullseye Urine Plate

Specialties

Urologists
General Practice
Pediatrics
OB/Gyn
Veterinarians

Cultures for

Enteric (E. coli, Klebsiella spp., Enterobacter spp., Proteus spp., Pseudomonas spp., Enterococci (Group D Streptococci), Staphylococci and Antibiotic Sensitivities



Southeast Vetlab Supply

18131 S.W. 98th Court
Miami, FL 33157-5509

305-253-1848 • 800-330-1522
Fax 305-232-8421